

# SENTINEL

3D PRINTED TITANIUM  
TUMOR CLIP MARKER

Available Codes	Measures	Q.ty per box
SNT 14/15	14G X CM.15	5
SNT 14/10	14G X CM.10	5
SNT 18/15	18G X CM.15	5
SNT 18/10	18G X CM.10	5
SNTM 14/15 (MRI)	14G X CM.15	5
SNTM 18/15 (MRI)	18G X CM.15	5



**TSUNAMI MEDICAL SRL**  
Via XXV Aprile n. 22  
41037 Mirandola (MO) Italy


UNI EN ISO 13485:2012

Tel. +39.0535.38397  
Fax. +39.0535.38399  
[www.tsunamimedical.it](http://www.tsunamimedical.it)  
[info@tsunamimed.com](mailto:info@tsunamimed.com)



# WHEN TO USE SENTINEL

- Before breast neoadjuvant therapy
- Before image guided radiotherapy of pancreas, liver and lung lesions
- During Brachytherapy for guiding conformal external beam
- After breast biopsy
- Before prostate radiotherapy to limit any damage in proximity to hollow organs



# FUTURES AND BENEFITS

- No migration thank's to the special shape and special 3d porous titanium surface
- 3D printed titanium realize a surface extremely compatible with ultrasound and easily visible
- Sentinel marker is suitable for MRI
- ergonomic handling with special soft touch polymers

